

Product datasheet for **HKRT502**

TissueScan, Kidney Cancer cDNA Array II

Product data:

Product Type:	TissueScan
Array Category:	Cancer Specific
Clinical Info:	/assets/documents/tissuescan/hkrt502.xls
Diagnosis Category:	Kidney Cancer
Tissue:	Kidney
Array Name:	Kidney Cancer cDNA Array II


Product images:



TissueScan product image.




[View online »](#)



A Human kidney TissueScan QPCR panel II was used in real-time PCR. The beta-actin contents in the 48 human cDNA were detected using SYBR Green I protocol and performed on ABI7900HT.

Fluorescence reading - cycle number



A Human kidney TissueScan QPCR panel II was used in SYBR Green I based real-time PCR. The expression profile of human Carbonic Anhydrase IX (CA9) gene (NM_001216) was detected using SYBR Green protocol and performed on ABI7900HT.

Normalized relative expression among the 48 kidney samples