

## Product datasheet for **HP200011**

### beta 2 Adrenergic Receptor (ADRB2) Human qPCR Primer Pair (NM\_000024)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	154
Forward Sequence:	TACCAGAGCCTGCTGACCAAGA
Reverse Sequence:	AGTCACAGCAGGTCTCATTGGC
ACCN:	<a href="#">NM_000024</a> , <a href="#">NM_000024.1</a> , <a href="#">NM_000024.2</a> , <a href="#">NM_000024.3</a> , <a href="#">NM_000024.4</a> , <a href="#">NM_000024.5</a> , <a href="#">BC012481</a> , <a href="#">BC012481.2</a> , <a href="#">BC063486</a> , <a href="#">BC073856</a> , <a href="#">NM_000024.6</a>
UniProt ID:	<a href="#">P07550</a>
Synonyms:	ADRB2R; ADRBR; B2AR; BAR; BETA2AR
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)