

Product datasheet for **HP200029**

Apolipoprotein H (APOH) Human qPCR Primer Pair (NM_000042)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	350
Forward Sequence:	AGAGTATGTCCTTTTCTGGAATC
Reverse Sequence:	GGCTCCATTTTCCTCCTCAGTG
ACCN:	NM_000042 , NM_000042.1 , NM_000042.2 , BC026283 , BC020703 , BG566874 , BM005736 , NM_000042.3
UniProt ID:	P02749
Synonyms:	B2G1; B2GP1; BG
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)