

## Product datasheet for **HP200030**

### CD95 (FAS) Human qPCR Primer Pair (NM\_000043)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	355
Forward Sequence:	GGACCCAGAATACCAAGTGCAG
Reverse Sequence:	GTTGCTGGTGAGTGTGCATTCC
ACCN:	<a href="#">NM_000043</a> , <a href="#">NM_000043.1</a> , <a href="#">NM_000043.2</a> , <a href="#">NM_000043.3</a> , <a href="#">NM_000043.4</a> , <a href="#">NM_000043.5</a> , <a href="#">BC012479</a> , <a href="#">BC012479.1</a> , <a href="#">BC065736</a> , <a href="#">BU620333</a>
UniProt ID:	<a href="#">P25445</a>
Synonyms:	ALPS1A; APO-1; APT1; CD95; FAS1; FASTM; TNFRSF6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)