

Product datasheet for **HP200031**

Androgen Receptor (AR) Human qPCR Primer Pair (NM_000044)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	367
Forward Sequence:	ATGGTGAGCAGAGTGCCCTATC
Reverse Sequence:	ATGGTCCCTGGCAGTCTCCAAA
ACCN:	<u>NM_000044</u> , <u>NM_000044.1</u> , <u>NM_000044.2</u> , <u>NM_000044.3</u> , <u>NM_000044.4</u> , <u>BC132975</u> , <u>NM_000044.6</u>
UniProt ID:	<u>P10275</u>
Synonyms:	AIS; AR8; DHTR; HUMARA; HYSY1; KD; NR3C4; SBMA; SMAX1; TFM
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)