

## Product datasheet for **HP200032**

### Arginase 1 (ARG1) Human qPCR Primer Pair (NM\_000045)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	383
Forward Sequence:	TCATCTGGGTGGATGCTCACAC
Reverse Sequence:	GAGAATCCTGGCACATCGGGAA
ACCN:	<a href="#">NM_000045</a> , <a href="#">NM_000045.1</a> , <a href="#">NM_000045.2</a> , <a href="#">NM_000045.3</a> , <a href="#">BC020653</a> , <a href="#">BC020653.1</a> , <a href="#">BC005321</a> , <a href="#">BG217880</a> , <a href="#">BG542163</a> , <a href="#">BT006741</a> , <a href="#">NM_000045.4</a>
UniProt ID:	<a href="#">P05089</a>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)