

## Product datasheet for **HP200033**

### ARSB Human qPCR Primer Pair (NM\_000046)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	411
Forward Sequence:	CACCAATGGCACAAAGCCTCTG
Reverse Sequence:	CGGTGAAGAGTCCACGAAGTTC
ACCN:	<a href="#">NM_000046</a> , <a href="#">NM_000046.1</a> , <a href="#">NM_000046.2</a> , <a href="#">NM_000046.3</a> , <a href="#">BC029051</a> , <a href="#">BC029051.1</a> , <a href="#">BC036609</a> , <a href="#">BM664932</a> , <a href="#">BQ216635</a> , <a href="#">BQ775217</a> , <a href="#">BX350102</a> , <a href="#">BX364854</a> , <a href="#">BX369127</a> , <a href="#">NM_000046.5</a>
UniProt ID:	<a href="#">P15848</a>
Synonyms:	ASB; G4S; MPS6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)