

Product datasheet for **HP200038**

ATM Human qPCR Primer Pair (NM_000051)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	472
Forward Sequence:	TGTTCCAGGACACGAAGGGAGA
Reverse Sequence:	CAGGGTTCTCAGCACTATGGGA
ACCN:	BC007023 , NM_000051 , NM_000051.3 , BC022306 , BC022307 , BC061584 , BC137169 , BC152385 , BC152389 , BT006764 , NM_000051.4
UniProt ID:	Q13315
Synonyms:	AT1; ATA; ATC; ATD; ATDC; ATE; TEL1; TELO1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)