

## Product datasheet for **HP200046**

### **BTD Human qPCR Primer Pair (NM\_000060)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	686
<b>Forward Sequence:</b>	TGGAGCCTCACCGCTTCAATGA
<b>Reverse Sequence:</b>	GGTCACTGCTATGACAAGGCTC
<b>ACCN:</b>	<a href="#">NM_000060</a> , <a href="#">NM_000060.1</a> , <a href="#">NM_000060.2</a> , <a href="#">NM_000060.3</a> , <a href="#">NM_000060.4</a> , <a href="#">BC012099</a> , <a href="#">BE814688</a> , <a href="#">BE814713</a> , <a href="#">BF195658</a> , <a href="#">BM565546</a> , <a href="#">BM710057</a> , <a href="#">BM839350</a> , <a href="#">BM919775</a> , <a href="#">BP326502</a>
<b>UniProt ID:</b>	<a href="#">P43251</a>
<b>Synonyms:</b>	biotinidase
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)