

Product datasheet for **HP200060**

CD40L (CD40LG) Human qPCR Primer Pair (NM_000074)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	959
Forward Sequence:	GCGGCACATGTCATAAGTGAGG
Reverse Sequence:	GTCCTTGCTTTTAACGGTCAGC
ACCN:	<u>NM_000074</u> , <u>NM_000074.1</u> , <u>NM_000074.2</u> , <u>BC071754</u> , <u>BC071754.1</u> , <u>BC074950</u> , <u>NM_000074.3</u>
UniProt ID:	<u>P29965</u>
Synonyms:	CD40L; CD154; gp39; hCD40L; HIGM1; IGM; IMD3; T-BAM; TNFSF5; TRAP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)