

Product datasheet for **HP200061**

CDK4 Human qPCR Primer Pair (NM_000075)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1019
Forward Sequence:	CCATCAGCACAGTTCGTGAGGT
Reverse Sequence:	TCAGTTCGGGATGTGGCACAGA
ACCN:	NM_000075 , NM_000075.1 , NM_000075.2 , NM_000075.3 , BC005864 , BC005864.2 , BC003644 , BC010153 , BC015669 , BI832658 , BQ773564 , NM_000075.4
UniProt ID:	P11802
Synonyms:	CMM3; PSK-J3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)