

## Product datasheet for **HP200090**

### CYP1B1 Human qPCR Primer Pair (NM\_000104)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1545
Forward Sequence:	GCCACTATCACTGACATCTTCGG
Reverse Sequence:	CACGACCTGATCCAATTCTGCC
ACCN:	<u><a href="#">NM_000104</a></u> , <u><a href="#">NM_000104.1</a></u> , <u><a href="#">NM_000104.2</a></u> , <u><a href="#">NM_000104.3</a></u> , <u><a href="#">BC012049</a></u> , <u><a href="#">BC012049.1</a></u> , <u><a href="#">BP326939</a></u> , <u><a href="#">BT019979</a></u> , <u><a href="#">BT020001</a></u>
UniProt ID:	<u><a href="#">Q16678</a></u>
Synonyms:	ASGD6; CP1B; CYP1B1; GLC3A; P4501B1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)