

Product datasheet for **HP200093**

Lipoamide Dehydrogenase (DLD) Human qPCR Primer Pair (NM_000108)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1738
Forward Sequence:	GGGACTAGAAGAGCTGGGAATTG
Reverse Sequence:	CATCCTCTGCTTTGTGAGCCAG
ACCN:	<u>NM_000108</u> , <u>NM_000108.1</u> , <u>NM_000108.2</u> , <u>NM_000108.3</u> , <u>NM_000108.4</u> , <u>BC018648</u> , <u>BC018648.2</u> , <u>BC018696</u> , <u>BI562331</u> , <u>NM_000108.5</u>
UniProt ID:	<u>P09622</u>
Synonyms:	DLDD; DLDH; E3; GCSL; LAD; OGDC-E3; PHE3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)