

Product datasheet for **HP200097**

SLC26A2 Human qPCR Primer Pair (NM_000112)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1836
Forward Sequence:	GCACTGTAACTTTATAGCTGGAG
Reverse Sequence:	AGTGAAGGAGGCACCACTGACA
ACCN:	<u>NM_000112</u> , <u>NM_000112.1</u> , <u>NM_000112.2</u> , <u>NM_000112.3</u> , <u>BC059390</u> , <u>BC059390.1</u> , <u>BI912879</u> , <u>BQ437711</u> , <u>BX640696</u>
UniProt ID:	<u>P50443</u>
Synonyms:	D5S1708; DTD; DTDST; EDM4; MST153; MSTP157
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)