

Product datasheet for **HP200106**

EPOR Human qPCR Primer Pair (NM_000121)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2057
Forward Sequence:	TTCTGTGCTTCACCGAGCGGTT
Reverse Sequence:	CGACACAGCTTCCATGGCTCAT
ACCN:	NM_000121 , NM_000121.1 , NM_000121.2 , NM_000121.3 , BC112153 , BC112153.1 , BC019092 , BC143697 , BC143698 , BC143700 , BU619338 , NM_000121.4
UniProt ID:	P19235
Synonyms:	EPO-R
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)