

Product datasheet for **HP200108**

XPG (ERCC5) Human qPCR Primer Pair (NM_000123)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2073
Forward Sequence:	GATGCGGAAGATTCGCTCCATG
Reverse Sequence:	GCTTTCCAGGAACATCTGTCCG
ACCN:	NM_000123 , NM_000123.1 , NM_000123.2 , NM_000123.3 , BC031522 , BC031522.2 , BE223050 , BQ002437 , BQ707436 , BX647399 , NM_000123.4
UniProt ID:	P28715
Synonyms:	COFS3; ERCC5-201; ERCM2; UVDR; XPG; XPGC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)