

Product datasheet for **HP200111**

ETFA Human qPCR Primer Pair (NM_000126)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 2108 |
| Forward Sequence: | TGGCGGTAGTGCCAGTTCAGAA |
| Reverse Sequence: | CTCTCCACTCTTCAAGCCTCGA |
| ACCN: | NM_000126 , NM_000126.1 , NM_000126.2 , NM_000126.3 , BC015526 , BC029995 , BC034156 , BC095457 , BG261160 , BT009796 , BU176141 , NM_000126.4 |
| UniProt ID: | P13804 |
| Synonyms: | EMA; GA2; MADD |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)