

Product datasheet for **HP200131**

Ferritin Light Chain (FTL) Human qPCR Primer Pair (NM_000146)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2512
Forward Sequence:	TACGAGCGTCTCCTGAAGATGC
Reverse Sequence:	GGTTCAGCTTTTCTCCAGGGC
ACCN:	<u>NM_000146</u> , <u>NM_000146.1</u> , <u>NM_000146.2</u> , <u>NM_000146.3</u> , <u>BC008439</u> , <u>BC008439.1</u> , <u>BC002991</u> , <u>BC004245</u> , <u>BC008441</u> , <u>BC013928</u> , <u>BC016346</u> , <u>BC016354</u> , <u>BC016715</u> , <u>BC018990</u> , <u>BC021670</u> , <u>BC058820</u> , <u>BC062708</u> , <u>BC067772</u> , <u>BX571748</u> , <u>NM_000146.4</u>
UniProt ID:	<u>P02792</u>
Synonyms:	LFTD; NBIA3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)