

## Product datasheet for **HP200153**

### GLI3 Human qPCR Primer Pair (NM\_000168)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2737
Forward Sequence:	TCAGCAAGTGGCTCCTATGGTC
Reverse Sequence:	GCTCTGTTGTCGGCTTAGGATC
ACCN:	<u><a href="#">NM_000168</a></u> , <u><a href="#">NM_000168.1</a></u> , <u><a href="#">NM_000168.2</a></u> , <u><a href="#">NM_000168.3</a></u> , <u><a href="#">NM_000168.4</a></u> , <u><a href="#">NM_000168.5</a></u> , <u><a href="#">BC117168</a></u> , <u><a href="#">BC032660</a></u> , <u><a href="#">BC113616</a></u> , <u><a href="#">BF194991</a></u> , <u><a href="#">BI497308</a></u> , <u><a href="#">NM_000168.6</a></u>
UniProt ID:	<u><a href="#">P10071</a></u>
Synonyms:	ACLS; GCPS; GLI3-190; GLI3FL; PAP-A; PAPA; PAPA1; PAPB; PHS; PPDIV
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)