

Product datasheet for **HP200158**

CD42b (GP1BA) Human qPCR Primer Pair (NM_000173)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2811
Forward Sequence:	ACCATCCTGGTGTCTGCCACAA
Reverse Sequence:	ACGGAGCTTTGGTGGCTGATCA
ACCN:	<u>NM_000173</u> , <u>NM_000173.1</u> , <u>NM_000173.2</u> , <u>NM_000173.3</u> , <u>NM_000173.4</u> , <u>NM_000173.5</u> , <u>NM_000173.6</u> , <u>BC027955</u> , <u>BE463855</u>
UniProt ID:	<u>P07359</u>
Synonyms:	BDPLT1; BDPLT3; BSS; CD42B; CD42b-alpha; DBPLT3; GP1B; GPIbA; GPIbalpha; VWDP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)