

Product datasheet for **HP200165**

GUCY2D Human qPCR Primer Pair (NM_000180)

Product data:

| | |
|--------------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 3000 |
| Forward Sequence: | GCCTCGTCTTTCTTGCTTCCT |
| Reverse Sequence: | GGAGAAAGGTGATGTCGTCCAC |
| ACCN: | <u>NM_000180</u> , <u>NM_000180.1</u> , <u>NM_000180.2</u> , <u>NM_000180.3</u> , <u>BC148421</u> , <u>BM931704</u> |
| UniProt ID: | <u>Q02846</u> |
| Synonyms: | CACD1; CG-E; CORD5; CORD6; CSNB1I; CYGD; GUC1A4; GUC2D; LCA; LCA1; RCD2; retGC; RETGC-1; ROS-GC1; ROSGC |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)