

## Product datasheet for **HP200171**

### Factor H (CFH) Human qPCR Primer Pair (NM\_000186)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3075
Forward Sequence:	ATTCCGTTCACTGCTACCACTTTG
Reverse Sequence:	CCATTGAGGAGTTCAGGAGGTG
ACCN:	<a href="#">NM_000186</a> , <a href="#">NM_000186.1</a> , <a href="#">NM_000186.2</a> , <a href="#">NM_000186.3</a> , <a href="#">BC012610</a> , <a href="#">BC037285</a> , <a href="#">BC073982</a> , <a href="#">BC110643</a> , <a href="#">BC142699</a> , <a href="#">BM842566</a> , <a href="#">BP324193</a> , <a href="#">NM_000186.4</a>
UniProt ID:	<a href="#">P08603</a>
Synonyms:	AHUS1; AMBP1; ARMD4; ARMS1; CFHL3; FH; FHL1; HF; HF1; HF2; HUS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)