

## Product datasheet for **HP200174**

### Hexokinase II (HK2) Human qPCR Primer Pair (NM\_000189)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3099
Forward Sequence:	GAGTTTGACCTGGATGTGGTTGC
Reverse Sequence:	CCTCCATGTAGCAGGCATTGCT
ACCN:	<a href="#">NM_000189</a> , <a href="#">NM_000189.1</a> , <a href="#">NM_000189.2</a> , <a href="#">NM_000189.3</a> , <a href="#">NM_000189.4</a> , <a href="#">BC064369</a> , <a href="#">BC064369.1</a> , <a href="#">BC021116</a> , <a href="#">BM706373</a> , <a href="#">BM912287</a> , <a href="#">NM_000189.5</a>
UniProt ID:	<a href="#">P52789</a>
Synonyms:	HKII; HXK2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)