

Product datasheet for **HP200197**

Integrin beta 4 (ITGB4) Human qPCR Primer Pair (NM_000213)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3691
Forward Sequence:	AGGATGACGACGAGAAGCAGCT
Reverse Sequence:	ACCGAGAACTCAGGCTGCTCAA
ACCN:	<u>NM_000213</u> , <u>NM_000213.1</u> , <u>NM_000213.2</u> , <u>NM_000213.3</u> , <u>NM_000213.4</u> , <u>BC109238</u> , <u>BC118916</u> , <u>BC126411</u> , <u>BC143738</u> , <u>BC143742</u> , <u>NM_000213.5</u>
UniProt ID:	<u>P16144</u>
Synonyms:	CD104; GP150
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)