

Product datasheet for **HP200239**

CMH1 (MYH7) Human qPCR Primer Pair (NM_000257)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4625
Forward Sequence:	GGAGTTCACACGCCTCAAAGAG
Reverse Sequence:	TCCTCAGCATCTGCCAGGTGT
ACCN:	<u>NM_000257</u> , <u>NM_000257.1</u> , <u>NM_000257.2</u> , <u>NM_000257.3</u> , <u>BC112171</u> , <u>BC112173</u> , <u>BF834726</u> , <u>NM_000257.4</u>
UniProt ID:	<u>P12883</u>
Synonyms:	CMD1S; CMH1; MPD1; MYHCB; SPMD; SPMM
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)