

## Product datasheet for **HP200256**

### ornithine aminotransferase (OAT) Human qPCR Primer Pair (NM\_000274)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4942
Forward Sequence:	CGTTGTCTGCTATCTCCAGTTCC
Reverse Sequence:	CCATGAACGCAGCCACATTTGG
ACCN:	<u>NM_000274, NM_000274.1, NM_000274.2, NM_000274.3, BC000964, BC016928, NM_000274.4</u>
UniProt ID:	<u>P04181</u>
Synonyms:	GACR; HOGA; OATASE; OKT
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)