

Product datasheet for **HP200258**

INPP5F (OCRL) Human qPCR Primer Pair (NM_000276)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 4952 |
| Forward Sequence: | TGGCACATCCTTAGAGGCTCTG |
| Reverse Sequence: | CACTGGCACCTTCATCCAAAGG |
| ACCN: | <u>NM_000276</u> , <u>NM_000276.1</u> , <u>NM_000276.2</u> , <u>NM_000276.3</u> , <u>BC094726</u> , <u>BC094726.1</u> , <u>BC018003</u> , <u>BC025253</u> , <u>BC130612</u> , <u>BC144106</u> , <u>NM_000276.4</u> |
| UniProt ID: | <u>Q01968</u> |
| Synonyms: | Dent-2; DENT2; INPP5F; LOCR; NPHL2; OCRL-1; OCRL1 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM. |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)