

Product datasheet for **HP200262**

PCBD1 Human qPCR Primer Pair (NM_000281)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5092
Forward Sequence:	CCGTGATGCCATCTTCAAGCAG
Reverse Sequence:	ATGGTGGTCCAGTTTCTCAGCC
ACCN:	<u>NM_000281</u> , <u>NM_000281.1</u> , <u>NM_000281.2</u> , <u>NM_000281.3</u> , <u>BC006324</u> , <u>BC006324.2</u> , <u>BI548743</u> , <u>BM724638</u> , <u>BQ685076</u> , <u>BQ897536</u> , <u>NM_000281.4</u>
UniProt ID:	<u>P61457</u>
Synonyms:	DCOH; PCBD; PCD; PHS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)