

Product datasheet for **HP200265**

PDHA1 Human qPCR Primer Pair (NM_000284)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5160
Forward Sequence:	GGATGGTGAACAGCAATCTTGCC
Reverse Sequence:	TCGCTGGAGTAGATGTGGTAGC
ACCN:	NM_000284 , NM_000284.1 , NM_000284.2 , NM_000284.3 , BC002406 , BC002406.1 , BU170389 , NM_000284.4
UniProt ID:	P08559
Synonyms:	PDHA; PDHAD; PDHCE1A; PHE1A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)