

Product datasheet for **HP200269**

PEX7 Human qPCR Primer Pair (NM_000288)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5191
Forward Sequence:	TCCTCATCACCTGTAGTGGCGA
Reverse Sequence:	GCTGTTACCTCTGGTTGGCT
ACCN:	NM_000288 , NM_000288.1 , NM_000288.2 , BC006268 , BC006268.1 , BC031606 , BP246238 , BQ011494 , NM_000288.4
UniProt ID:	O00628
Synonyms:	PBD9B; PTS2R; RCDP1; RD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)