

Product datasheet for **HP200281**

Phospholipase A2 IIA (PLA2G2A) Human qPCR Primer Pair (NM_000300)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5320
Forward Sequence:	GTTGACGACAGGAAAGGAAGCC
Reverse Sequence:	CCAGACGTTTGTAGCAACAGTCA
ACCN:	<u>NM_000300</u> , <u>NM_000300.1</u> , <u>NM_000300.2</u> , <u>NM_000300.3</u> , <u>BC005919</u> , <u>BM708748</u> , <u>BU626291</u> , <u>BX641714</u> , <u>NM_000300.4</u>
UniProt ID:	<u>P14555</u>
Synonyms:	MOM1; PLA2; PLA2B; PLA2L; PLA2S; PLAS1; sPLA2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)