

Product datasheet for **HP200306**

PITX2 Human qPCR Primer Pair (NM_000325)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5308
Forward Sequence:	GTGTGGACCAACCTTACGGAAG
Reverse Sequence:	CGAAGCCATTCTTGCATAGCTCG
ACCN:	<u>NM_000325</u> , <u>NM_000325.1</u> , <u>NM_000325.3</u> , <u>NM_000325.5</u> , <u>BC013998</u> , <u>BC013998.2</u> , <u>BC106010</u> , <u>BC106010.1</u> , <u>BP372081</u> , <u>BX380370</u> , <u>NM_000325.6</u>
UniProt ID:	<u>Q99697</u>
Synonyms:	ARP1; ASGD4; Brx1; IDG2; IGDS; IGDS2; IHG2; IRID2; Otlx2; PTX2; RGS; RIEG; RIEG1; RS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)