

## Product datasheet for **HP200318**

### delta Sarcoglycan (SGCD) Human qPCR Primer Pair (NM\_000337)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6444
Forward Sequence:	GACAGAGCTGAGACTGGAATCC
Reverse Sequence:	GAAGACCTTCTGCCTCGTTCCT
ACCN:	<a href="#">NM_000337</a> , <a href="#">NM_000337.1</a> , <a href="#">NM_000337.2</a> , <a href="#">NM_000337.3</a> , <a href="#">NM_000337.4</a> , <a href="#">NM_000337.5</a> , <a href="#">BC020740</a> , <a href="#">BC020740.1</a> , <a href="#">BC008431</a> , <a href="#">BQ434817</a> , <a href="#">BQ446019</a> , <a href="#">BX501169</a> , <a href="#">BX537948</a> , <a href="#">BX642431</a>
UniProt ID:	<a href="#">Q92629</a>
Synonyms:	35DAG; CMD1L; DAGD; LGMDR6; SG-delta; SGCDP; SGD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)