

Product datasheet for **HP200322**

SLC3A1 Human qPCR Primer Pair (NM_000341)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 6519 |
| Forward Sequence: | TACGGTTCTGGCTCACAAAGGG |
| Reverse Sequence: | GCTCCGAGTATTGTGTGACCGT |
| ACCN: | <u>NM_000341</u> , <u>NM_000341.1</u> , <u>NM_000341.2</u> , <u>NM_000341.3</u> , <u>BC022386</u> , <u>BC022386.1</u> , <u>BC093624</u> , <u>BC093626</u> |
| UniProt ID: | <u>Q07837</u> |
| Synonyms: | ATR1; CSNU1; D2H; NBAT; RBAT |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)