

## Product datasheet for **HP200327**

### SOX9 Human qPCR Primer Pair (NM\_000346)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6662
Forward Sequence:	AGGAAGCTCGCGGACCACTAC
Reverse Sequence:	GGTGGTCCTTCTGTGCTGCAC
ACCN:	<u><a href="#">NM_000346</a></u> , <u><a href="#">NM_000346.1</a></u> , <u><a href="#">NM_000346.2</a></u> , <u><a href="#">NM_000346.3</a></u> , <u><a href="#">BC056420</a></u> , <u><a href="#">BC056420.1</a></u> , <u><a href="#">BC007951</a></u> , <u><a href="#">BC018276</a></u> , <u><a href="#">BT006875</a></u> , <u><a href="#">NM_000346.4</a></u>
UniProt ID:	<u><a href="#">P48436</a></u>
Synonyms:	CMD1; CMPD1; SRA1; SRXX2; SRXY10
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)