

## Product datasheet for **HP200333**

### **SUR1 (ABCC8) Human qPCR Primer Pair (NM\_000352)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	6833
<b>Forward Sequence:</b>	GACGACAAGAGGACAGTGGTCT
<b>Reverse Sequence:</b>	GCATTCAGACCTCTGGAAGTCC
<b>ACCN:</b>	<u><a href="#">NM_000352</a></u> , <u><a href="#">NM_000352.1</a></u> , <u><a href="#">NM_000352.2</a></u> , <u><a href="#">NM_000352.3</a></u> , <u><a href="#">NM_000352.4</a></u> , <u><a href="#">BC167829</a></u> , <u><a href="#">BM313666</a></u> , <u><a href="#">BM352492</a></u> , <u><a href="#">BU950000</a></u> , <u><a href="#">NM_000352.6</a></u>
<b>UniProt ID:</b>	<u><a href="#">Q09428</a></u>
<b>Synonyms:</b>	ABC36; HHF1; HI; HRINS; MRP8; PHHI; PNDM3; SUR; SUR1; SUR1delta2; TNDM2
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)