

## Product datasheet for **HP200338**

### TGFBI Human qPCR Primer Pair (NM\_000358)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7045
Forward Sequence:	GGACATGCTCACTATCAACGGG
Reverse Sequence:	CTGTGGACACATCAGACTCTGC
ACCN:	<u><a href="#">NM_000358</a></u> , <u><a href="#">NM_000358.1</a></u> , <u><a href="#">NM_000358.2</a></u> , <u><a href="#">BC000097</a></u> , <u><a href="#">BC000097.2</a></u> , <u><a href="#">BC004972</a></u> , <u><a href="#">BC026352</a></u> , <u><a href="#">BT009820</a></u> , <u><a href="#">NM_000358.3</a></u>
UniProt ID:	<u><a href="#">Q15582</a></u>
Synonyms:	BIGH3; CDB1; CDG2; CDGG1; CSD; CSD1; CSD2; CSD3; EBMD; LCD1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)