

## Product datasheet for **HP200355**

### UROS Human qPCR Primer Pair (NM\_000375)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7390
Forward Sequence:	CAGACAGTTGCACACCCAGGAA
Reverse Sequence:	GTGCTTGAGACTGTATGTGAGGC
ACCN:	<u><a href="#">NM_000375</a></u> , <u><a href="#">NM_000375.1</a></u> , <u><a href="#">NM_000375.2</a></u> , <u><a href="#">BC002573</a></u> , <u><a href="#">BC002573.2</a></u> , <u><a href="#">BC004338</a></u> , <u><a href="#">BC024936</a></u> , <u><a href="#">BI668388</a></u> , <u><a href="#">NM_000375.3</a></u>
UniProt ID:	<u><a href="#">P10746</a></u>
Synonyms:	UROIIIS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)