

Product datasheet for **HP200356**

Vitamin D Receptor (VDR) Human qPCR Primer Pair (NM_000376)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7421
Forward Sequence:	CGCATCATTGCCATACTGCTGG
Reverse Sequence:	CCACCATCATTACACGAACTGG
ACCN:	NM_000376 , NM_000376.1 , NM_000376.2 , BC060832 , BC060832.1 , BC033465 , BM685276 , BM908433 , BQ002808 , BX355584 , NM_000376.3
UniProt ID:	P11473
Synonyms:	NR111; PPP1R163
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)