

Product datasheet for **HP200399**

CD41 (ITGA2B) Human qPCR Primer Pair (NM_000419)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3674
Forward Sequence:	CTGTCCAGCTACTGGTGCAAGA
Reverse Sequence:	ATGTTGTGCCCCAGTGGCTCCAA
ACCN:	<u>NM_000419</u> , <u>NM_000419.1</u> , <u>NM_000419.2</u> , <u>NM_000419.3</u> , <u>NM_000419.4</u> , <u>BC117443</u> , <u>BC126442</u> , <u>BE409167</u> , <u>BM148153</u> , <u>NM_000419.5</u>
UniProt ID:	<u>P08514</u>
Synonyms:	BDPLT2; BDPLT16; CD41; CD41B; GP2B; GPIIb; GT; GT1; GTA; HPA3; PPP1R93
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)