

Product datasheet for **HP200403**

Cytokeratin 2 (KRT2) Human qPCR Primer Pair (NM_000423)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3849
Forward Sequence:	ACCTACCGCAAAGCTGGAGG
Reverse Sequence:	CAGAACCTCCAAAGGCAGCCTT
ACCN:	NM_000423 , NM_000423.1 , NM_000423.2 , BC096294 , BC099643 , BC099644 , NM_000423.3
UniProt ID:	P35908
Synonyms:	CK-2e; K2e; KRT2A; KRT2E; KRTE
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)