

Product datasheet for **HP200405**

L1CAM Human qPCR Primer Pair (NM_000425)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3897
Forward Sequence:	TCGCCCTATGTCCACTACACCT
Reverse Sequence:	ATCCACAGGGTTCTTCTCTGGG
ACCN:	NM_000425 , NM_000425.1 , NM_000425.2 , NM_000425.3 , NM_000425.4 , BC126229 , BC025843 , BC136447 , NM_000425.5
UniProt ID:	P32004
Synonyms:	CAML1; CD171; HSAS; HSAS1; MASA; MIC5; N-CAM-L1; N-CAML1; NCAM-L1; S10; SPG1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)