

Product datasheet for **HP200412**

Myosin Light Chain 2 (MYL2) Human qPCR Primer Pair (NM_000432)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4633
Forward Sequence:	CGGAGAGGTTTTCCAAGGAGGA
Reverse Sequence:	CTCTTCTCCGTGGGTGATGATG
ACCN:	BC015821 , NM_000432 , NM_000432.1 , NM_000432.2 , NM_000432.3 , BC015821.1 , BC031006 , BC031008 , BX538165 , NM_000432.4
UniProt ID:	P10916
Synonyms:	CMH10; MLC-2s/v; MLC2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)