

Product datasheet for **HP200420**

PDE6 alpha (PDE6A) Human qPCR Primer Pair (NM_000440)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5145
Forward Sequence:	GTCCGTGCTTTCCTCAACTGTG
Reverse Sequence:	GGACCAGAGTAAGGTGGAAGCTTC
ACCN:	NM_000440 , NM_000440.1 , NM_000440.2 , BC035909 , BC035909.1 , BC144044 , BU734142 , BX537592 , BX648986 , NM_000440.3
UniProt ID:	P16499
Synonyms:	CGPR-A; PDEA; RP43
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)