

## Product datasheet for **HP200422**

### CD31 (PECAM1) Human qPCR Primer Pair (NM\_000442)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5175
Forward Sequence:	AAGTGGAGTCCAGCCGCATATC
Reverse Sequence:	ATGGAGCAGGACAGGTTCAATC
ACCN:	<u><a href="#">NM_000442</a></u> , <u><a href="#">NM_000442.1</a></u> , <u><a href="#">NM_000442.3</a></u> , <u><a href="#">NM_000442.4</a></u> , <u><a href="#">BC051822</a></u> , <u><a href="#">BC051822.1</a></u> , <u><a href="#">BC022512</a></u> , <u><a href="#">BC062608</a></u> , <u><a href="#">BG272085</a></u>
UniProt ID:	<u><a href="#">P16284</a></u>
Synonyms:	CD31; CD31/EndoCAM; endoCAM; GPIIA'; PECA1; PECAM-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)