

Product datasheet for **HP200425**

Plectin (PLEC) Human qPCR Primer Pair (NM_000445)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5339
Forward Sequence:	AGCGTGAGAAGGAGAAGCTCCA
Reverse Sequence:	CAGAGAGGAAGCTTTGCTGCAG
ACCN:	<u>NM_000445</u> , <u>NM_000445.1</u> , <u>NM_000445.2</u> , <u>NM_000445.3</u> , <u>NM_000445.4</u> , <u>BC007597</u> , <u>BC013206</u> , <u>BM681271</u> , <u>BP231761</u> , <u>BP312608</u> , <u>BP315878</u> , <u>BP321327</u> , <u>BP359323</u> , <u>NM_000445.5</u>
UniProt ID:	<u>Q15149</u>
Synonyms:	EBS1; EBSMD; EBSND; EBSO; EBSOG; EBSPA; HD1; LGMD2Q; LGMDR17; PCN; PLEC1; PLEC1b; PLTN
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)