

## Product datasheet for **HP200438**

### **HNF1 beta (HNF1B) Human qPCR Primer Pair (NM\_000458)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	6928
<b>Forward Sequence:</b>	CCCAGCAAATCTTGACCAGGC
<b>Reverse Sequence:</b>	ACCTCAGTGACCAAGTTGGAGC
<b>ACCN:</b>	<u><a href="#">NM_000458</a></u> , <u><a href="#">NM_000458.1</a></u> , <u><a href="#">NM_000458.2</a></u> , <u><a href="#">NM_000458.3</a></u> , <u><a href="#">BC017714</a></u> , <u><a href="#">BC017714.1</a></u> , <u><a href="#">BT007126</a></u> , <u><a href="#">NM_000458.4</a></u>
<b>UniProt ID:</b>	<u><a href="#">P35680</a></u>
<b>Synonyms:</b>	ADTKD3; FJHN; HNF-1-beta; HNF-1B; HNF1beta; HNF2; HPC11; LF-B3; LFB3; MODY5; RCAD; T2D; TCF-2; TCF2; VHNF1
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)