

Product datasheet for **HP200440**

Thrombopoietin (THPO) Human qPCR Primer Pair (NM_000460)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7066
Forward Sequence:	CCAGAGGTTCACCCCTTGCCTA
Reverse Sequence:	CCAGAATGTCCTGTGCCTTGGT
ACCN:	NM_000460 , NM_000460.1 , NM_000460.2 , NM_000460.3 , BC130322 , BC136338 , BC143982 , BX108191 , NM_000460.4
UniProt ID:	P40225
Synonyms:	MGDF; MKCSF; ML; MPLLG; THCYT1; TPO
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)