

Product datasheet for **HP200446**

Twist (TWIST1) Human qPCR Primer Pair (NM_000474)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7291
Forward Sequence:	GCCAGGTACATCGACTTCCTCT
Reverse Sequence:	TCCATCCTCCAGACCGAGAAGG
ACCN:	<u>NM_000474</u> , <u>NM_000474.1</u> , <u>NM_000474.2</u> , <u>NM_000474.3</u> , <u>BC036704</u> , <u>BC036704.2</u> , <u>NM_000474.4</u>
UniProt ID:	<u>Q15672</u>
Synonyms:	ACS3; bHLHa38; BPES2; BPES3; CRS; CRS1; CSO; SCS; SWCOS; TWIST
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)