

Product datasheet for **HP200494**

HOXA13 Human qPCR Primer Pair (NM_000522)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3209
Forward Sequence:	TGGAACGGCCAAATGTACTGCC
Reverse Sequence:	GGTATAAGGCACGCGCTTCTTTC
ACCN:	NM_000522 , NM_000522.1 , NM_000522.2 , NM_000522.3 , NM_000522.4 , BC075791 , BC075791.1 , BC016712 , NM_000522.5
UniProt ID:	P31271
Synonyms:	HOX1; HOX1J
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)